

CNS Vital Signs Advancing Multiple Sclerosis Care

**Adding Value to Your Practice by Providing Solutions for Measuring,
Monitoring and Managing Neurocognitive and Behavioral Health...**

CNS
Vital Signs®

www.CNSVS.com

Contents

■ Why CNS Vital Signs?	3
■ Why CNS Vital Signs in Multiple Sclerosis?.....	4
■ About CNS Vital Signs Assessment Platforms	6
■ CNS Vital Signs Multiple Sclerosis Toolbox	9
■ Optimized for MS Evaluation & Management.....	11
■ Helping Assess Comorbidities	13
■ HOW? Brief Interpretation	16
■ How...Next Steps	19

The following pages have been assembled from various sources and publications and is meant to be a reference or roadmap guide to assist and inform how CNS Vital Signs can be used to improve clinical insight and care management, enable current guidelines, be integrated into a clinic or practice, and help improved practice revenues and performance.



Why CNS Vital Signs?


CNS Vital Signs **valid, reliable, and affordable** 'research quality' NEUROCOGNITIVE & BEHAVIORAL HEALTH assessment platform can be easily configured and deployed depending on each practices or researchers needs and goals. The CNS Vital Signs assessment platforms helps to support a practices comprehensive, state-of-the-art clinical assessment, and evidence-based treatment services for children, adolescents, and adults across the lifespan by:

- Accurately measuring and characterizing a patient's neurocognitive function based on his or her status or effort
- Facilitating the thinking about the patient's condition (50+ well known medical and health rating scales) and helping to explain the patient's current difficulties
- Optimizing serial administration which helps to monitor and guide effective intervention.
- *Systematically collecting brain function, behavioral, symptom, and comorbidity data enabling outcomes and evidence-based medicine*



Enhanced Brain & Behavior Evaluation and Care Management

OBJECTIVE, PRECISE, and STANDARDIZED... Customizable Toolboxes or Test Panels Supporting many Neurological, Psychiatric, & Psychological Clinical Guidelines



Extend Practice Efficiency

Objective and Evidence-Based Assessments, Auto-Scored and Systematically Documented.

(HIPAA Enabled)

Enhanced Revenue Streams

Expanded Services with Well Established Billing Codes to **Improve Practice Referrals** and Performance



Why CNS Vital Signs in Multiple Sclerosis?

Benefits for Multiple Sclerosis Clinics

National Multiple Sclerosis Society: Expert Opinion Paper Summary

■ Assessment and Management of Cognitive Impairment in Multiple Sclerosis

- Cognitive deficits appear to be present in more than half of MS patients, *however the majority of persons with MS do not have impairments that significantly impair daily functioning*
- Learning/memory, speed of information processing, working memory, cognitive flexibility and other executive functions appear to be most commonly impaired
- ***Periodic screening for such deficits is recommended.***
- Intervention for such deficits is recommended: Training in strategies to compensate for deficits, Counseling / psychotherapy for patients and family to address accompanying behavioral changes and emotional responses, and develop realistic expectations
- Treatment with medications (disease-modifying and/or symptomatic therapies)

■ Enhanced MS Evaluation, Management & Tracking Strategies

- ***CNS Vital Signs provides a valid, reliable and granular view of neurocognitive status***
- **Efficient:** Reports are Auto-Scored in seconds and Screens for possible in-valid tests
- **Multi-Modal Assessment platform** allowing for improved **Comorbid Symptom identification and management** e.g. Fatigue, Depression, Mood, Quality of Life / Outcomes, Etc.
- **Longitudinal reports** auto-generated to **monitor and measure** e.g. treatment outcomes

■ Increased Revenues

The standard for specific assessment of cognitive function in MS patients has been the comprehensive neuropsychological assessment. Adding CNS Vital Signs for the early detection, characterization, and monitoring of MS cognitive dysfunction progression should be part of routine care as an in-take baseline, as part of a full neuropsychological assessment and periodic retesting providing clinicians with a valid and reliable longitudinal view that can be beneficial both clinically and in counseling patients and working with family members.



HOW Can CNS VS Help?

Where the CNS Vital Signs Assessment Platform can add efficiency, validity, and standardization to the assessment of brain functions, symptoms, and treatment effects.

Etiology of MS

Genes & Environment


Abnormal Immune Response

Multiple Sclerosis

Immune Stimulation (Infections?)

Secondary Factors

Worsening MS



Adapted From: Cohen, B. Int J MS Care, 2007; 9 (Suppl4):4-10; Ben-Zacharia A, et al. Multiple Sclerosis Counseling Points. 2007;3(3):1-12.

Adapted From: Crayton H, et.al. Neurology, 2004;63 (11 Suppl 5): S12-S18.



About CNS Vital Signs?

Assessing Brain Function: CNS Vital Signs is a clinical testing procedure used by clinicians to evaluate and manage the neurocognitive state of a patient. Across the lifetime, serial testing allows ongoing assessments of a patient's condition, disease progression, or clinical outcome.

About CNS Vital Signs

Both Valid & Reliable Neurocognitive Testing and Evidence-Based Symptom & Functional Ratings Scales in one Platform

Optimized for...

- **MULTI-MODAL Assessment** enabling the *efficient collection* and *systematic documentation* of important brain function and behavioral, symptom and comorbid clinical endpoints
- **Lifespan Testing** - *Rapid Neurocognitive Testing from ages 8 to 90*
- **Longitudinal View** - *CNS Vital Signs contains an Auto-Randomization Algorithm... Ideal for Serial Neurocognitive Testing with an almost unlimited number of alternate forms (others use a pseudo-randomization or limited number of alternate forms)*
- **Flexible Deployment** - *Easy Integration via Local Computer Software and Web-Based Testing Solutions... Ideal for busy clinics, hospitals, or academic research*

Clinician Benefits

- **RAPID INSIGHT...** computerized neurocognitive testing helps clinicians evaluate and describe the health of the cognitive or higher functions of the brain in a more granular and standardized fashion.
- **DASHBOARD VIEW...** Neurocognitive domain functions and functional status is presented in a summary view that is easy to interpret.
- **LONGITUDINAL VIEW...** Repeated testing allows clinicians to track disease progress and treatment/rehabilitation effects
- **DETAILED VIEW...** Each report presents the testing data in a detailed view. All results can be easily exported to EMR's or spreadsheets for clinical or research purposes.
- **VALID ACROSS the LIFE SPAN...** Peer reviewed normative data allows clinicians to examine patients from age 8 to 90.



About: CNS Vital Signs in Multiple Sclerosis

CNS Vital Signs provides clinicians and researchers with leading edge neurocognitive and behavioral health assessment technologies that efficiently collect valid and reliable brain & behavioral clinical endpoints for a more objective view of a patient's functional status, disease progression, and outcomes. The CNS Vital Signs Assessment platform supports a lifespan chronic care model and helps enable productive interactions between the family, caregivers, and a specialist practice team.

CNS Vital Signs is a clinical procedure that utilizes scientifically validated objective tests to evaluate the neurocognitive status of patients and covers a range of mental processes from simple motor performance, attention and memory, to executive functions. The CNS Vital Signs tests are computerized versions of well established neuropsychological tests. Medical professionals and researchers know that good health has many dimensions, one of the most important and yet least measured is the health of a person's brain. Outcomes based medicine seeks a quantitative estimate of the effect of impairment or disease and the effectiveness and efficiency of treatment. CNS Vital Signs provides a **standardized** and **quantitative** view of your patient's **CORE COGNITIVE FUNCTION**.

CNS Vital Signs computerized neuropsychological tests can enhance efficiency and insight in assessing cognitive status and the difference between "normal performance" and a patient's current status and provides the clinician with a normative comparison that can be paired with an interview, exam, and other valid test(s) or rating scales to help add validity to the evaluation and management of Multiple Sclerosis. Re-evaluation or serial testing with CNS Vital Signs supports the effective management and tailoring of treatments e.g., medications and assessment of outcomes. A very detailed assessment of abilities is auto-scored, and the pattern of strengths and weaknesses can be used in treatment planning and measuring progress.

One of the most robust features of the CNS Vital Signs assessment is its randomization algorithm allowing for an almost infinite number of alternate forms. This allows for retesting patients and minimal practice effects. Clinicians establish a baseline and upon re-test, compare the results to assist in decision-making regarding the observed change in the patient's condition, monitor disease or recovery progress, measure treatment results, compliance, and outcomes e.g., Therapy Management, Medication Optimization, Etc. Often Patients and families benefit from seeing testing results allowing the understanding of the status and nature of their or a loved one's neurocognitive function. CNS Vital Signs is one of many tools clinicians use in evaluating changes in a patient's condition.

If you have a question or would like to register for a free in-service webinar go to www.CNSVS.com or email support@cnsvs.com or call 1.888.750.6941.



Why Use CNS Vital Signs to Assess MS?

The CNS Vital Signs VSX Assessment Platform represents a legacy of innovation and a commitment to advancing neurocognitive and behavioral clinical assessment tools.

Clinical Pathology

Measure and Monitor

Assess **BRAIN FUNCTION** and Determine the Existence or Level of **IMPAIRMENT...**

CNS Vital Signs computerized neurocognitive testing allows clinicians to **assess abnormal neurocognitive impairment** by comparing patients to a 'PEER REVIEWED' normative data set from **ages 8 to 90** across the lifespan

Provides a broad spectrum of clinical domains and the sensitivity to assess neurocognitive function to reveal abnormal cognitive function.

Comorbid Status

Measure and Monitor

Assess symptoms or **COMORBID conditions...**

Evidence-based rating scales and neurocognitive testing can help clinicians **sort out symptom, behavioral, and comorbid issues** and help better understand possible brain and behavior relationships.

50+ Free Rating Scales:

- SF – 36 Medical Outcomes
- Zung Self-Rating Anxiety and Depression Scales
- NeuroPsych Questionnaire NPQ-207 & NPQ-45 both Child & Adult

Serial Assessment

Longitudinal View

KEY ADVANTAGE

...contains an **auto-randomization algorithm...** Ideal for serial testing with an **almost unlimited number of alternate forms** (other systems use a pseudo-randomization or limited number of alternate forms).

This allows practices to shift toward new assessment approaches that allow for monitoring of change and the reinforcement of treatment compliance.

Neurocognitive function in MS patients can vary and affect various sites within the brain, no single cognitive deficit pattern is characteristic of all patients with MS. However, specific deficits tend to be seen in MS. These deficits include problems with memory, attention, executive functioning, generative fluency, and information processing speed. Adapted from: Rao et.al. ; Heaton et.al.; Benedict et.al.



Advancing MS Care Management

CNS Vital Signs Multiple Sclerosis Toolbox



Clinician Expertise

Brain Function: Processing Speed,
Memory, Attentional, Executive,
Psychomotor Speed & more

Behaviors,
Symptoms, and
Comorbidities

Computerized Neurocognitive Testing

- Nine Neurocognitive Domains Measured
- Processing and Psychomotor Speed
- Frontal Lobe / Executive & Attentional Tests
- Recognition Memory – immediate and delayed recall
- Immediate Auto – Scored Reports
- Rapid Assessment – 30 -45 Minute initial Assessment/Baseline, 15 - 45 Minute for monitoring
- Easy to interpret and longitudinally graph
- Systematic & Standardized Documentation for Patient Registry/Research
- HIPAA Compliant

Computerized Medical and Health Rating Scales*


- SF – 36 Medical Outcomes
- Zung Self-Rating Anxiety and Depression Scales
- NeuroPsych Questionnaire NPQ-207 & NPQ-45 both Child & Adult
- Epworth Sleepiness
- Pittsburgh Sleep Quality Index

* Used with permission... Free use of rating scales




CNS Vital Signs MS Toolbox

MOS SF-36... Widely Used Measure



Optimized for MS Evaluation & Management



CNS Vital Signs Neurocognitive Battery in Multiple Sclerosis

Results: Utilizing data from all 42 patients together, there was a diffuse pattern of cognitive impairment compared to age-matched controls in all cognitive domains tested ($p < 0.02$). However, when divided into high and low functioning groups, the high functioning group had a more specific cognitive pattern, with particular difficulties with complex information processing (symbol digit coding, shifting attention test) and working memory. The low functioning group continued to have a diffuse impairment pattern.

Conclusions/Relevance: With a cognitively high functioning group of RRMS patients with well controlled MS, a subcortical pattern emerges, with particular difficulties with complex information processing and working memory. The cognitive pattern is much more diffuse with the low functioning group, even after controlling for motor speed and overall reaction time. These results could help explain the variance in cognitive testing that can be seen in MS patients.

Adapted from: AAN 2009; Higher Cognitively Functioning Relapsing-Remitting Multiple Sclerosis Patients Have a More Specific Pattern of Impairment on Neuropsychological Testing Sandeep Vaishnavi, MD, PhD, John Barkenbus, MD, C. Thomas Gualtieri, MD; NC Neuropsychiatry; Raleigh & Charlotte, NC

The CNS Vital Signs Sleep Toolbox helps clinicians systematically collect brain function, symptoms, and comorbidities data, automatically scoring and systematically documenting the resulting clinical endpoints.



Computerized neuropsychological battery detects psychomotor processing speed impairment as the core cognitive deficit in relapsing - remitting multiple sclerosis patients

Messinis L., PhD., Anyfantis E., M.Sc., Lyros E., MD, PhD., Papathanasopoulos P., MD, PhD Department of Neurology, Neuropsychology Section, University of Patras Medical School, Greece; 7th International Congress on the Improvement of the Quality of Life on Dementia, etc; Thessaloniki, Greece, 29 January 1st February 2009

Discussion and Conclusions

Our findings support the literature that RRMS patients present cognitive impairments when compared to healthy controls [7]. Furthermore, measures of information processing speed appear to be the most robust and sensitive markers of this impairment [8], a finding confirmed by our data.

These results also indicate that the CNS vital signs battery is able to detect cognitive impairments in RRMS patients. Furthermore, it appears to be highly sensitive in detecting psychomotor processing speed deficits in RRMS patients with a low disability status, as were the patients in this study.

We conclude, that RRMS patients present cognitive deficits even in the absence of significant physical disability, and that these patients should be routinely assessed for cognitive decline. The CNS Vital signs neuropsychological battery could provide the non - neuropsychologist clinician with a reliable screening tool for detecting cognitive deficits in RRMS patients.

Table 2. Neuropsychological performance on the CNSVS Battery: mean (SD)

Cognitive Domain	RRMS	Controls	
Memory	93.8 (8.5)	94.7 (5.8)	NS
*Psychomotor Speed	155.23 (20.68)	177.86 (22.02)	< .05
Reaction Time	682.06 (115.08)	658.97 (82.25)	NS
*Complex Attention	8.8 (6.2)	5.8 (2.6)	< .05
*Cognitive Flexibility	38.4 (15.2)	45.6 (10.10)	< .05

*Significant difference between groups (p <.05) based on parametric t-tests

Table 1. clinical and demographic characteristics: mean (SD)

Variables	RRMS (n= 30)	Controls (n=25)
*Age	34.6 (8.4)	35.4 (7.6)
*Sex	27 females / 3 males	23 females/ 2 males
* Education (yrs)	12.2 (3.0)	12.7 (2.6)
* § Estimated Intelligence level	98.50 (12.30)	99.40 (11.65)
Disease duration	5.8 (4.6)	-
EDSS (disability status)	2.8 (1.4)	-

*Non significant difference between groups (p <.05) based on parametric t-tests (age, education, intelligence level) and χ^2 test (sex)


§ Intelligence (IQ) was estimated by administering the vocabulary and matrix reasoning subscales of the Wechsler abbreviated scale of intelligence (WASI), Greek-adapted version [6]



Tools to Help Assess Symptoms and Comorbidities

NPQ – 45

Rapid In-take or Re-test to Assess the Neuro–Psych Status of a Patient



The *Neuropsych Questionnaire (NPQ) Short Form (SF - 45)* provides a subjective measure of 13 neuropsych symptoms. The symptoms are **Attention, Impulsive, Memory, Anxiety, Panic, Depression, Mood Stability, Oppositional (child – adolescent), Aggression, Fatigue, Sleep, Suicide, and Pain**. The shorter NPQ version is used to monitor or follow-up with the patient before or during their visit. The NPQ 45 can be used when the longer version is either impractical or inappropriate e.g. the physician wants a quick view of their patients core symptoms. **Both versions are automatically scored and the data stored.**

NPQ – 45 (Adult Patient & Informant Version)


Rapid In-take or Re-test to Assess the Neuro–Psych Status of a Patient

NeuroPsych Questionnaire (NPQ) SF-45 (Page 1 of 2)			
Subject Reference/ID: NPQ45SymptomAdult		Test Date: March 29, 2009 15:35:40	
Age: 64		Administrator: Neuropsych Solutions	
Total Test Time: 0:28 (min:secs) for all tests in this report		Language: English (United States)	
<i>This scale was administered using CNS Vital Signs</i>			
Domain	Score	Severity	Description
Attention	100	Mild	The Neuropsych Questionnaire Short Form asks patients (or an appropriate observer) a series of questions about their clinical state. The questions are about the symptoms of various neuropsychiatric disorders. The terminology is similar to that used in the diagnostic manuals, and in many familiar clinical questionnaires and rating scales; but it has been simplified, and all symptoms are scored on the same metric. Scores are reported on a scale of 0 (not a problem) to 300 (severe). As a rule, scores above 225 indicate a severe problem; scores from 150-224 indicate a moderate problem; and scores from 75-149, a mild problem. A high score on the Neuropsych Questionnaire Short Form means that the patient is reporting more symptoms of greater intensity. It doesn't necessarily mean that the patient has a particular condition; just that he or she (or their spouse, parent or caregiver) are saying that they have a lot of intense symptoms. Conversely, a low score simply means that the patient (or caregiver) is not reporting symptoms associated with a particular condition, at least during the period of time specified. It does not mean that the patient does not have the condition. Just as some people over-state their problems, others tend to under-state their problems. The Neuropsych Questionnaire Short Form is not a diagnostic instrument. The results it generates are only meant to be interpreted by an experienced clinician in the course of a clinical examination.
Impulsive	160	Moderate	
Memory	125	Mild	
Anxiety	167	Moderate	
Panic	100	Mild	
Depression	160	Moderate	
Mood Stability	125	Mild	
Aggression	200	Moderate	
Fatigue	167	Moderate	
Sleep	100	Mild	
Suicide	100	Mild	
Pain	225	Severe	
Attention Questions			
1	Difficulty concentrating	1 - A mild problem	
2	Easily distracted	1 - A mild problem	
3	Feeling scattered, disorganized	1 - A mild problem	
4	Forgetful, I need constant reminding	0 - Not a problem	
5	Short attention span	2 - A moderate problem	
Impulsive Questions			
1	Feeling restless	3 - A severe problem	
2	Fidgety, I can't sit still	1 - A mild problem	
3	Impatient	3 - A severe problem	
4	Impulsive, act without thinking	0 - Not a problem	
5	Overly active	1 - A mild problem	
Memory Questions			
1	Forgetful, I need constant reminding	0 - Not a problem	
2	My mind goes blank	2 - A moderate problem	
3	Problems with memory	3 - A severe problem	
4	Putting something down and then forgetting where you put it	0 - Not a problem	
Anxiety Questions			
1	Feeling anxious	0 - Not a problem	
2	Feeling nervous	2 - A moderate problem	
3	Feeling restless	3 - A severe problem	
4	Feeling tense	2 - A moderate problem	
5	Fidgety, I can't sit still	1 - A mild problem	
6	Worrying too much	2 - A moderate problem	
Panic Questions			
1	Attacks of intense anxiety	0 - Not a problem	
2	Feeling so nervous it's hard to breathe	2 - A moderate problem	
3	Panic attacks	1 - A mild problem	
Depression Questions			
1	Feeling depressed	1 - A mild problem	
2	Feeling discouraged about the future	3 - A severe problem	
3	Feeling irritable	0 - Not a problem	
4	Feeling little or no interest in things	1 - A mild problem	
5	Not enjoying things as much as before	3 - A severe problem	
Mood Stability Questions			
1	Anger	3 - A severe problem	
2	Easily frustrated	2 - A moderate problem	
3	Feeling irritable	0 - Not a problem	
4	My moods change quickly	0 - Not a problem	


Key MS Symptoms



Computerized Neurocognitive Testing can Help Assess Cognitive Function vs. Possible Comorbidities



Cognition and Depression



“Indeed, there is some suggestion that cognitive or executive functioning deficits may be a trait risk factor for depression (Douglas and Porter, 2009; Frasch et al., 2009; Micco et al., 2009; Reppermund et al., 2009). Furthermore, worse neuropsychological test performance at baseline is associated with poorer response to treatment (Dunkin et al., 2000; Kampf-Sherf et al., 2004; Mohlman and Gorman, 2005), and cognitive deficits are more pronounced in patients who are unemployed (Baune et al., 2010). It is possible that treatment refractory depression is a subtype characterized in part by cognitive impairment.

The accurate identification and quantification of neurocognitive impairment are important for research relating to neurobiological underpinnings, treatment, and functional outcome in patients with mood disorders. It is essential, methodologically, that we have accurate methods for identifying those patients who are objectively cognitively impaired and separate them from patients who have the subjective experience of poor thinking skills or thinking that is easily perturbed by negative affect, but perform normally on cognitive testing in controlled conditions. The treatments and outcomes for these two groups may differ markedly, as well as the prognosis.”

Source: Identifying a cognitive impairment subgroup in adults with mood disorders. J Affect Disord. 2011 Aug;132(3):360-7. Epub 2011 Mar 25.

<http://www.ncbi.nlm.nih.gov/pubmed/21439647>



Cognition and Depression

Cognitive Flexibility

Domain scored from two venerable AD/HD tests

- Healthy Control
- Mood Disorder, Normal Cognition
- ◆ *Mood Disorder, Cognitive Impairment*





Fig. 3. Distributions of CNS Vital Signs *cognitive flexibility* index score in patients with or without impaired cognition. Figure note: Healthy control, N=660. Mood disorder, normal cognition, n=128. Mood disorder, cognitive impairment, n=58. *Normative scores were truncated at 40. Each value represents the percentage of subjects in that score range.



HOW? CNS Vital Signs begins with...



A: Conducting a Valid Assessment (Refer to the Test Administration Guide.) To begin the staff should collect information about the CHIEF or REFERRAL COMPLAINT. This will be a primary driver for the selection of tests and rating scales. For initial evaluations or in complex presentations, a broad spectrum battery is always an appropriate starting point.

B: Review the immediately auto-scored report to **1** validate testing effort, **2** evaluate the Domain Dashboard to quickly assess the level of impairment or grade the deficit, and **3** Evaluate the Domain Pattern to help rule-in, rule-out, or confirm certain clinical conditions. Feedback to the patient on the testing results may be presented at the clinical encounter or at a subsequent patient visit.

C: If invalid test results were noted then consider re-testing the patient to confirm clinical results. If the test results were valid, then, as part a continuum of care, reschedule testing to track disease progression and measure ongoing status or outcomes.


NOTE: The **Validity Indicator** denotes a guideline for representing the possibility of an invalid test or domain score. “No” means a clinician should evaluate whether or not the test subject understood the test, put forth their best effort, or has a clinical condition requiring further evaluation.



Evaluate Severity – Impairment Status

CNS Vital Signs grades *severity of impairment* based on an age-matched normative comparison database... mTBI Example

Psychometric and Normative Comparison



Above:	> 110	High Function and High Capacity
Average:	90 - 110	Normal Function and Normal Capacity
Low Average:	80 - 90	Slight Deficit and Slight Impairment
Low:	70 - 79	Moderate Deficit and Impairment Possible
Very Low:	< 70	Deficit and Impairment Likely


Standard Scores

CNS Vital Signs Clinical Report					Test Date: July 23 2011 10:48:38				
Subject ID: mTBI or AD/HD					Administrator: Technician				
Language: English (United States)					Age: 27				
Patient Profile:	Percentile Range				> 74	25 - 74	9 - 24	2 - 8	< 2
	Standard Score Range				> 109	90 - 109	80 - 89	70 - 79	< 70
Domain Scores	Subject Score	Standard Score	Percentile	VI**	Above	Average	Low Average	Low	Very Low
Neurocognition Index (NCI)	NA	85	16	Yes			X		
Composite Memory	102	103	58	Yes		X			
Verbal Memory	51	93	32	Yes		X			
Visual Memory	51	110	75	Yes	X				
Processing Speed	48	79	8	Yes				X	
Executive Function	34	75	5	Yes				X	
Psychomotor Speed	174	93	32	Yes		X			
Reaction Time*	555	107	68	Yes		X			
Complex Attention*	21	56	1	Yes					X
Cognitive Flexibility	26	63	1	Yes					X
Total Test Time (min: secs)	29:12				Total time taken to complete the tests shown.				
Domain Dashboard: Above average domain scores indicate a standard score (SS) greater than 109 or a Percentile Rank (PR) greater than 74, indicating a high functioning test subject. Average is a SS 90-109 or PR 25-74, indicating normal function. Low Average is a SS 80-89 or PR 9-24 indicating a slight deficit or impairment. Below Average is a SS 70-79 or PR 2-8, indicating a moderate level of deficit or impairment. Very Low is a SS less than 70 or a PR less than 2, indicating a deficit and impairment. Reaction times are in milliseconds. An * denotes that "lower is better", otherwise higher scores are better. Subject Scores are raw scores calculations generated from data values of the individual subtests.									
VI** - Validity Indicator: Denotes a guideline for representing the possibility of an invalid test or domain score. "No" means a clinician should evaluate whether or not the test subject understood the test, put forth their best effort, or has a clinical condition requiring further evaluation.									
Verbal Memory Test (VBM)	Score	Standard	Percentile	Verbal Memory Test: Subjects have to remember 15 words and recognize them in a field of 15 distractors. The test is repeated at the end of the battery. The VBM test measures how well a subject can recognize, remember, and retrieve words e.g. exploit or attend literal representations or attribute. "Correct Hits" refers to the number of target words recognized. Low scores indicate verbal memory impairment.					
Correct Hits - Immediate	13	102	55						
Correct Passes - Immediate	14	95	37						
Correct Hits - Delay	9	85	16						
Correct Passes - Delay	15	109	73						
Visual Memory Test (VIM)	Score	Standard	Percentile	Visual Memory Test: Subjects have to remember 15 geometric figures, and recognize them in a field of 15 distractors. The test is repeated at the end of the battery. The VIM test measures how well a subject can recognize, remember, and retrieve geometric figures e.g. exploit or attend symbolic or spatial representations. "Correct Hits" refers to the number of target figures recognized. Low scores indicate visual memory impairment.					
Correct Hits - Immediate	13	107	68						
Correct Passes - Immediate	14	117	87						
Correct Hits - Delay	13	111	77						
Correct Passes - Delay	11	93	32						
Finger Tapping Test (FTT)	Score	Standard	Percentile	The FTT is a test of motor speed and fine motor control ability. There are three rounds of tapping with each hand. The FTT test measures the speed and the number of finger-taps with each hand. Low scores indicate motor slowing. Speed of manual motor activity varies with handedness. Most people are faster with their preferred hand but not always.					
Right Taps Average	64	104	61						
Left Taps Average	60	105	63						
Symbol Digit Coding (SDC)	Score	Standard	Percentile	The SDC test measures speed of processing and draw upon several cognitive processes simultaneously, such as visual scanning, visual perception, visual memory, and motor functions. Errors may be due to impulsive responding, misperception, or confusion.					
Correct Responses	50	80	9						
Errors*	2	92	30						
Stroop Test (ST)	Score	Standard	Percentile	The ST measures simple and complex reaction time, inhibition / disinhibition, mental flexibility or directed attention. The ST helps assess how well a subject is able to adapt to rapidly changing and increasingly complex set of directions. Prolonged reaction times indicate cognitive slowing / impairment. Errors may be due to impulsive responding, misperception, or confusion.					
Simple Reaction Time*	231	108	70						
Complex Reaction Time Correct*	542	100	50						
Stroop Reaction Time Correct*	568	112	79						
Stroop Commission Errors*	8	5	1						
Shifting Attention Test (SAT)	Score	Standard	Percentile	The SAT measures executive function or how well a subject recognizes set shifting (mental flexibility) and abstraction (rules, categories) and manages multiple tasks simultaneously. Subjects have to adjust their responses to randomly changing rules. The best scores are high correct responses, few errors and a short reaction time. Normal subjects may be slow but accurate, or fast but not so accurate. Attention deficit may be apparent.					
Correct Responses	47	82	12						
Errors*	13	75	5						
Correct Reaction Time*	1003	97	42						
Continuous Performance Test (CPT)	Score	Standard	Percentile	The CPT measures sustained attention or vigilance and choice reaction time. Most normal subjects obtain near-perfect scores on this test. A long response time may suggest cognitive slowing and/or impairment. More than 2 errors (total) may be clinically significant. More than 4 errors (total) indicate attentional dysfunction.					
Correct Responses	40	104	61						
Omission Errors*	0	104	61						
Commission Errors*	0	108	70						
Choice Reaction Time Correct*	400	99	47						



Neurocognitive Domain Dashboard

CNS Vital Signs presents testing results in Subject (raw), Standard Scores, and Percentile Ranks. NOTE: See the CNS Vital Signs Interpretation Guide for more information.



Patient Profile:	Percentile Range				> 74	25 - 74	9 - 24	2 - 8	< 2
	Standard Score Range				> 109	90 - 109	80 - 89	70 - 79	< 70
Domain Scores	Subject Score	Standard Score	Percentile	VI**	Above	Average	Low Average	Low	Very Low
Neurocognition Index (NCI)	NA	85	16	Yes			x		
Composite Memory	102	103	58	Yes		x			
Verbal Memory	51	93	32	Yes		x			
Visual Memory	51	110	75	Yes	x				
Processing Speed	48	79	8	Yes				x	
Executive Function	34	75	5	Yes				x	
Psychomotor Speed	174	93	32	Yes		x			
Reaction Time*	555	107	68	Yes		x			
Complex Attention*	21	56	1	Yes					x
Cognitive Flexibility	26	63	1	Yes					x
Total Test Time (min: secs)	29:12				Total time taken to complete the tests shown.				

SD = Standard Deviation from the MEAN



HOW can CNS Vital Signs Benefit My Practice?

Ask about our NO COST Practice Evaluation!

CNS Vital Signs Benefits



Enhanced Patient Insight and Care Management



Enables Evidence-Based Medicine and Outcomes



Improved Practice Efficiencies and Documentation



Improved Practice Revenues and Performance

Solution Example

CNS Vital Signs Mobile Test Station ULTRA Series



\$1,400.00
Testing Station
with 40 test
sessions.

Potential Return On Investment

Based on Established Billing Codes*

40 Patient Test Sessions ROI:

\$2,400 to \$10,000+

**Possible Yearly IMPACT... \$80K to \$160K
depending on patient volumes...**

*Based on a survey of Payers. Contact support@cnsvs.com for billing information.

Popular with Clinics and Hospitals: Engineered with BUSY PRACTICES in mind (roll into exam rooms), the Ultra Series combines the ultimate in practical functionality, ergonomic ease-of-use, and remarkable durability.



Solutions for Measuring, Monitoring, and
Managing Neurocognitive and Behavioral Health



NEXT STEPS:

Contact Us...

Getting Started

Step One: Register at www.CNSVS.com
After registering download the VSX 'Brief-Core'
Assessment Software with 5 FREE Test Sessions...
Take it for a test drive.

Step Two: *Schedule a FREE One-on-One In-Service Webinar...* Contact CNS Vital Signs Support support@cnsvs.com with dates and times that you will be available.

After the webinar the total CNS Vital Signs Assessment platform (Web & Local) can be configured to meet your practice needs.

Learn More

Contact me to receive report examples, case studies, administration guides etc.

- **Website:** www.CNSVS.com
- **Phone:** 888.750.6941
- **Email:** support@cnsvs.com
- **Address:**
598 Airport Blvd.
Suite 1400
Morrisville, NC 27560

"The webinar training was terrific... it covered the Validity & Reliability of the platform, the interpretation of results, billing and coding, testing protocol, and the integration of the CNS Vital Signs platform into our practice." Practice Administrator

